

5. CHAPTER V

PREPARATION OF SOLUTIONS

Use accurately calibrated equipment, which means meets the TS specifications.

Working temperature of standard solution should approximate that of its temperature during preparation and standardisation. If temperature correction necessary, sufficient accuracy may be obtained by use of Tables.

Table 5.1. Temperature corrections for volume of aqueous solutions (AOAC, 1990).

Volume of standard solution (ml)	Correction ml at						
	8 °C	12 °C	16 °C	20 °C	24 °C	28 °C	30 °C
10	0.01	0.01	0.01	0.00	0.00	-0.02	-0.02
20	0.03	0.02	0.01	0.00	-0.02	-0.03	-0.03
25	0.03	0.03	0.02	0.00	-0.02	-0.04	-0.05
30	0.04	0.03	0.02	0.00	-0.02	-0.05	-0.07
40	0.06	0.04	0.03	0.00	-0.03	-0.07	-0.09
50	0.07	0.06	0.03	0.00	-0.04	-0.09	-0.12

5.1. Preparations of Indicators

5.1.1. Phenolphthalein

1 % phenolphthalein (w/w) solution is prepared by using alcohol.

5.1.2. Methylene Blue

1 % methylene blue (w/w) solution is prepared by using water.

5.1.3. Methylene Orange

0.1 % methylene orange (w/w) solution is prepared by using water.

5.1.4. Ferric Alum

Saturated solution of $\text{FeNH}_4(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$ in water.

5.1.5. Potassium Chromate

5 % potassium chromate, K_2CrO_4 , (w/w) solution is prepared by using water.

5.1.6. Methylene Red

1 % methylene red solution is prepared by using water.

5.2. Standard Solution of Hydrochloric Acid, HCl

Table 5.2. Approximate volumes of 36.5–38 % HCl required to make 10 liter standard solution.

Approximate normality	ml HCl to be diluted to 10 liter
0.01	8.6
0.02	17.2
0.10	86.0
0.50	430.1
1.00	860.1

Standardisation of HCl Solutions

Titrate 40 ml against alkaline solution (NaOH) of same concentration as acid being standardised in 300 ml flask that has been swept free from CO_2 , using CO_2 -free water and 3 drops of phenolphthalein as indicator.

Normality = (ml standard alkali * normality of alkali) / ml HCl

If more concentrated than desired, dilute solution to required normality value by following formula:

$$V_1 = V_2 * N_2 / N_1$$

where N_2 and V_2 represent normality and volume of stock solution and V_1 volume of stock solution should be diluted to obtain desired normality, N_1 .

Check the normality of final solution by titration as above.

5.3. Preparation of Standard Sodium Hydroxide Solution, NaOH

To 1 part NaOH (reagent quality containing < 5 % Na_2CO_3) in flask add 1 part water and swirl until solution is complete.

Table 5.3 Approximate volumes of NaOH solution necessary to make 10 liter standard solution.

Approximate normality	ml NaOH to be diluted to 10 liter
0.01	8.6
0.02	17.2
0.10	86.0
0.50	430.1
1.00	860.1

Add required volume of NaOH solution (1+1) to 10 liter CO₂-free water. Check normality, as in standardisation section, which should be slightly high and adjust to desired concentration by following formula

$$V_1 = V_2 * N_2 / N_1$$

where N₂ and V₂ represent normality and volume of stock solution and V₁ volume of stock solution should be diluted to obtain desired normality, N₁.

Check the normality of final solution by titration as above.

Standardisation of NaOH Solution

Accurately weigh enough dried acid potassium phthalate, KHC₈H₄O₄, to titrate 40 ml and transfer to 300 ml flask that has been swept free from CO₂. Add 50 ml cool CO₂-free water. Stopper flask and swirl gently until sample dissolved. Titrate to pH 8.6 with solution being standardised, taking 3 drops of phenolphthalein as indicator. Determine volume NaOH required to produce end point of blank by matching color in another flask containing 3 drops of phenolphthalein and same volume of CO₂-free water. Subtract volume required from that used in first titration and calculate normality.

$$\text{Normality} = \frac{\text{g KHC}_8\text{H}_4\text{O}_4}{\text{ml NaOH}} * \frac{1000}{204.229}$$

5.4. Preparation of Standard Sulphuric Acid Solution, H₂SO₄

Table 5.4. Approximate volumes of 95–98 % H₂SO₄ required to make 10 liter standard solution.

Approximate normality	ml H ₂ SO ₄ to be diluted to 10 liter
0.01	2.8
0.02	5.6
0.10	27.7
0.50	138.1
1.00	276.1

Standardisation of H₂SO₄ Solutions

Standardisation of H₂SO₄ can be carried out as in HCl (Section 5.2)

5.5. Standard Solution of Ammonium and Potassium Thiocyanates

Prepare 0.1 N solution from reagents free from Cl; using 7.612 g NH₄SCN in 1 liter or 9.718 g KSCN in 1 liter.

5.6. Standard Solution of Sodium Thiosulphate

Dissolve 25 g NaS₂O₃·5H₂O in 1 liter of water. Boil gently 5 min. transfer while hot to storage bottle previously cleaned with hot chromic acid cleaning solution and rinsed with warm boiled water.

5.7. Standard Solution of Silver Nitrate

Dissolve slightly more than theoretical weigh of AgNO₃ (equivalent weight, 169.87) in halogen-free water, and dilute to volume.

Standardisation by Mohr Method:

Potassium chloride: Recrystallize KCl 3 times from water, dry at 110 °C, and then heat at 500 °C to constant weight. Equivalent weight KCl = 74.555.

Potassium chromate: 5 % solution of K₂CrO₄ in water.

Accurately weigh enough KCl to yield titration of 40 ml (0.3 g for 0.1 N solution) and transfer to 250 ml Erlenmeyer with 40 ml water. Add 1 ml K₂CrO₄ solution and titrate with AgNO₃ until first perceptible pale-red-brown appears. From titration volume subtract ml of the AgNO₃ solution required to produce end point color in 75 ml water containing 1 ml K₂CrO₄ solution. From the volume of AgNO₃ calculate normality:

$$\text{Normality} = \text{g KCl} * 1000 / \text{ml AgNO}_3 * 74.555$$

5.8. Standard Solution of Iodine

Dissolve weighed amounts of I (127.1 g/L) and KI, in proportion of 20 g KI to 13 g I in 50 ml water, when I dissolves, transfer solution to volumetric flask. Dilute to volume with water and mix thoroughly. Store in dark brown bottle in dark place.

5.9. Standard Solution of Alcoholic Potassium Hydroxide

Preparation of 0.5 N KOH: Dissolve 28.055 g KOH in 1 liter of alcohol.

5.10. Standard Solution of Sodium Carbonate

Preparation of 0.5 M NaCO₃: Dissolve 5.3 g NaCO₃ in water and dilute to 100 ml.

5.11. Standard Solution of Boric Acid

Preparation of 0.2 M H₃BO₃: Dry H₃BO₃ to constant weight in desiccator over CaCl₂. Dissolve 21.405 g H₃BO₃ in water and dilute to 1 liter.

5.12. Standard Solutions of Fehling's

I: Dissolve 69.278 g copper sulphate CuSO₄.5H₂O in water make up to 1 liter.

II: Dissolve 100 g NaOH and 346 g sodium potassium tartrate in water and make up to 1 liter.

5.13. Standard Starch Solution

Add desired powdered starch into volumetric flask. Add some water. Place the volumetric flask into water bath at 80 °C and swirl frequently until starch dissolved. Cool the flask to 20 °C and make up to desired volume.