

# FE 204 Experiment 1

## Preparation of Bacteriological Media and Sterilization

# Contents

- Media and types of media
- Preparation of General Purpose Media
- Procedure
- Sterilization techniques
- Working principle of autoclave

# Media and Types of Media

**Definition:** Any material which supports the growth and replication of microorganisms is called **Medium (pl. media)**.

## • Liquid Media

Generally called as **Broth**.  
For ex: Nutrient broth

## • Solid media

Generally called as **Agar**.  
For ex: Nutrient agar

Any broth media can be turned into a solid media by adding agar agar into broth with a 2.5% ratio of volume.

# Special Purpose Media

Anaerobic media

Synthetic media

Transport media

Enriched media

Selective media

Differential media

Selective & Differential media

Microbiological assay media

- **Anaerobic media:** Media supporting growth of only anaerobic m/os O<sub>2</sub> is removed from the media.
- **Synthetic media:** referred to as chemically defined medium which their concentrations are known.
- **Transport media:** are used for temporary storage maintaining the number of m/os during transportation.
- **Enriched media:** contain nutrients which supports all types of microorganisms.
- **Selective media:** support the growth of a certain type of m/o while inhibiting the growth of others.
- **Differential media:** contain indicators that distinguish between organisms on the basis of their appearance on the medium.
- **Selective & Differential media:** hold the properties of both selective and differential media.
- **Microbiological assay media:** are used to measure concentrations of substances such as antibiotics and vitamins. For ex: Muller-Hinton Agar.

## Sterilization Methods



Sterilization by Heat

Sterilization by Radiation

Sterilization by Filtration

Chemical Sterilization

# Autoclave

- For sterilization of materials and media
- 121 C for sterilization of media
- 134 C for sterilization of waste and dish



# Oven and Incubators

- Dry air sterilization Oven 180 C



- Incubator
- 25, 37 and 45 C





# Materials used in LAB

- 250 ml beaker
- 500 ml flask
- 10 ml pipettes
- Sterile petri dishes
- Nutrient agar
- Nutrient broth
- Sensitive balance
- Magnetic stirrer
- Peptone
- Electrical heater
- Clean spatula



# Preperation Materials

- Balance
- Spatules and spoons




# Pipette and Petri Dishes




# Magnetic Stirrer

- Mixing and heating solutions homogenously.



# Procedure



Medium	Amounts	Container	Sterilization
a) Broth	5-10 ml	Test Tube	In Autoclave
b) Agar Slants	10 ml	Test Tube	In Autoclave
c) Agar Deep	7 ml	Test Tube	In Autoclave
d) Agar Petri Plates	20 – 25 ml	Petri Dish	Media → In autoclave Petri Dishes → In Oven

# Procedure


- Preparation of Nutrient broth:
- Calculate ingredients to make 100 ml of Nutrient broth by looking the formula given on the box of media.

## Preparation

Suspend 25 g Standard I Nutrient Broth/litre, autoclave (15 min at 121 °C).

pH:  $7.5 \pm 0.2$  at 25 °C.

The prepared media are clear and yellowish-brown.





# Procedure

- Preparation of Nutrient broth:
- Weigh the required calculated amount on balance and put inside a 250 ml beaker
- Add 100 ml distilled water and mix by magnetic stirrer until getting a clear solution without any suspended powder.





# Procedure

- Preparation of Nutrient broth:
- Dispense into test tubes of 5ml amount and place inside a tube-rack.
- Sterilize all of the tubes in autoclave at 121 C for 20 min. After sterilization they are ready to use.



# Procedure

- Preparation of Nutrient agar:
- Calculate ingredients to make **200 ml of Nutrient agar** by looking the formula given on the box of media.
- Weigh the required calculated amount on balance and put inside a 250 ml balloon/flask.
- Add 200 ml of distilled water, place magnetic stirrer and mix the ingredients and heat by electrical heater.



Agar Slant



Agar Deep



Agar Petri Plate

Heat  
until  
boiling



# Procedure

- Preparation of Nutrient agar:
- Continue to mix and shake until boiling the whole suspension in the flask. Careful about not spilling around the neck.
- After a clear solution obtained. Distribute into test tubes for;
- Agar slant and agar deep
- Keep rest in flask for petri plates.
- Sterilize in autoclave at 121 C for 20 min.



Mix  
and  
heat



Agar Slant

Pour 10 ml  
of hot  
molten  
agar for  
agar slant



Agar Deep

Pour 7 ml  
of hot  
molten  
agar for  
agar deep.


# Procedure

After sterilization of agar tubes and flask.

- For Agar slant;
  - Keep the tubes at slant position lay down horizontally at room temperature to solidify.
- For Agar Deep;
  - Keep the tubes at vertical position until solidify.
- For Agar Petri Plates;
  - Pour 20-25 ml of Molten Agar in Sterile Petri Plates at aseptic conditions and keep them at room temperature to solidify.



# Result



Medium	Amounts	Container	Sterilization
a) Broth	5-10 ml	Test Tube	In Autoclave
b) Agar Slants	10 ml	Test Tube	In Autoclave
c) Agar Deep	7 ml	Test Tube	In Autoclave
d) Agar Petri Plates	20 – 25 ml	Petri Dish	Media → In autoclave Petri Dishes → In Oven