# FE 305 Experiment 2 Bacteriological Examination of Water

### Contents

- Total Aerobic Mesophilic Bacteria Count
- Coliform Count
  - Presumptive Test
  - Confirmed Test
  - Completed Test

#### Samples needed for this experiment

- River Water directly obtained
- Pool water
- Tap water
- Bottled Drinking water

Experiment No from Book: 4-15-46

### Water types

- Spring or well water
- Specifically prepared drinking water
- Purified water
- Fluoridated water
- Mineral water

### Importance of examination of water

- Main food material
- Can be consumed directly or added to food materials
- In order to keep product safe added water should be safe



# Preperation of Serial Dilution

- In order to make a succesfull counting method we need serial dilution of samples.
- 9 ml of Sterile 0.1 % peptone water tubes are used for dilution.
- 1 ml sample was taken from non diluted sample and added into 9 ml of dilution water in order to obtain 10<sup>-1</sup> dilution.



# Total Aerobic Mesophilic Bacteria Count

• Media: Plate Count Agar



- Incubation conditions:
- Spread Plate Technique
- Count formed colonies



### Material and Methods

- Sterile PCA petri plates
- Sterile Serial Dilution water tubes (9 ml)
- Spreader
- Alcohol
- Pipette (5 and 1 ml)
- Bunsen burner
- Incubator 37 C

# Procedure of Total Count

- Take 0.1 ml of sample from each dilution by using sterile pipette at Aseptic Conditions.
- Place the sample on PCA petri plate near flame.





# Procedure of Total Count Spread Plate

- Sterilize spreader with dipping in alcohol and passing through flame.
- Spread the sample on petri plate by using spreader.



 Incubate inoculated petri dishes at 37°C Oven for 24-48h. Count formed colonies

# Results of Total Count

• Count formed colonies on every petri plate for each dilution and record results on the given table below.



- You can divide petri into equal parts to count easily.
- Ex: For river water
  - Non  $\rightarrow$  Too nuber to count  $\rightarrow$  TNTC • -1  $\rightarrow$  546
  - -2 → 186
  - -3 → 34
  - -4 → 0

Sample names	Dilutions						
	non	10-1	10-2	10 <sup>-3</sup>	10-4	10 <sup>-5</sup>	
River water	TNTC	546	186	34	0	0	
Pool water							
Tap water							
Drinking water							

# Calculation

- # of microorganisms / g or ml of water= $\frac{Count of microorganisms in one petri}{inoculum amount in one petri} X Dilution factor$
- Ex: For river water take only counts between 30<x<300 colonies into calculation.
  - Non: TNTC 546  $\rightarrow$  No calculation above 300 • -1: 186 186 • -2: 0.1 *ml* 34 34 • -3: 0.1 ml $\rightarrow$  No calculation below 30 0 • -4:
  - Take average of these two results
  - # of microorganisms / g or ml of water=  $\frac{186000+340000}{2}$  = 263000 microorganism / ml of water

# Coliform Count

- Presumptive Test
  - Media: Brillant Green Bile Broth
  - Incubation: 37 C for 48 hours
  - Most Probable Number method 3 tubes
  - Decide Positive tubes and calculate # of coliforms / ml water
- Confirmed Test
  - Media: Eosine Methylene Blue Agar
  - From Positive Presumptive Test Tubes
  - Streak Plate Technique
  - Observe isolated colonies on petri
- Completed Test
  - Samples from colonies formed on petri plates
  - Observation under microscope 40x and 100x with Gram Staining

# Presumptive Test MPN method materials

- Water sample
- Serial Dilution Water Test tubes
- BGBB with durham tubes
- Sterile pipette (10 and 1 ml)
- Bunsen burner
- Tube rack

# MPN method procedure

- Take 1ml of sample from each dilution by sterile pipette and inoculate 1ml in each of three seperate BGBB broth tubes at aseptic conditions.
- By this, you will have 3 BGBB tubes for each dilution sample inoculated.



- Place the tubes in a tube rack in order from non diluted to diluted.
- Incubate at 37°C for 48 h. Observe results and decide positive test tubes.



### Results of MPN method



Gas formation in the durham tubes

Turbidity at the broth



2 1 0 from MPN Table

- For each dilution count positive test tubes from 3 test tubes.
- Record results for each dilution.
- Calculate # of coliforms / ml of water from MPN table

Sample	non	<b>10</b> <sup>-1</sup>	<b>10</b> <sup>-2</sup>	<b>10</b> -3	10-4	<b>10</b> <sup>-5</sup>
River water	3	3	2	1	0	0

### Calculation

		Confidence limits (95%)		# of coliforms / ml of water =	
Positive tubes	MPN/g ou ml	Low	High	MPN number from table	
0-0-0	<3.0		9.5	100 × Dilution Factor of middle tube	
0-0-1	3.0	0.15	9.6		
0-1-0	3.0	0.15	11		
0-1-1	6.1	1.2	18	• For 2 1 0 $\rightarrow$ MPN number from table = 15	
0-2-0	6.2	1.2	18		
0-3-0	9.4	3.6	38		
1-0-0	3.6	0.17	18	15	
1-0-1	7.2	1.3	18	$\frac{1}{100}$ X Dilution factor of middle tube	
1-0-2	11	3.6	38	100	
1-1-0	7.4	1.3	20		
1-1-1	11	3.6	38		
1-2-0	11	3.6	42	$\frown$	
1-2-1	15	4.5	42	2 (1) 0	
1-3-0	16	4.5	42		
2-0-0	9.2	1.4	38	<ul> <li>Look at the records at dilution table and find for dilution of</li> </ul>	
2-0-1	14	3.6	42	10-3	
2-0-2	20 MIDN	4.5	42	middle tube which is $10^{-3} \rightarrow \text{DIIUTION Factor } 10^{-3}$	
2-1-0	15	3.7	42	1 🗗	
2-1-1	20 number	4.5	42	$\frac{15}{10}$ X 103	
2-1-2	27	8.7	94	100 1 105	
				150 coliforms / ml of water	

### Confirmed test

- Made from positive presumptive test tubes
- To confirm culture growing in tubes as coliforms

#### **Materials**

- Sterile EMBA plates
- Positive presumptive test tubes
- Inoculating loops
- Bunsen burner

#### **Method**

• Streak Plate Method

#### Incubation

• Incubate inoculated petri plates at 37°C for 24 – 48 hours.

### Procedure

- Sterilize the inoculating loops by holdingon flame of bunsen burner.
- Take 1-2 loops of culture sample from Positive tubes
- Inoculate sample by using streak plate method and incubate.

#### The Streak Plate Isolation Method



### Results





### Completed test

- Morphological examination aand Gram staining of growing cultures on EMBA.
- Most of the lactose fermenters seems as Gram negative (-)

#### Materials

- Light Microscope with high power and oil immersion objectives
- Gram staining kit (Crystal Violet, Safranine, Iodine, Alcohol solutions)
- Pasteur pipette
- Glass slide
- Forceps
- Inoculating loop

#### Method

• Gram Staining

### Procedure

 Take samples from EMBA plate by using sterile loop and place on a sterile glass slide.



# **Smear Preparation**



Allow the smear to air dry



Pass through flame several times to kill and fix microorganisms to the slide

# Staining Smear



### Observation Under Microscope

 After staining place glass slides on Microscope stage and put a drop of immersion oil. Observe under 100x objective and decide G + or G – according to colors.





# Evaluation of Results

• Total count → Calculate and record # of microorganisms / ml of water

#### Coliform count

- Presumptive test MPN method
  - Calculate and record # of coliforms / ml of water (It should be zero for drinking water
- Confirmed test
  - Record EMBA growth results of positive presumptive test tubes.
- Completed test
  - Record Gram stain results of colonies isolated on EMBA plate

Presence of coliforms in drinking water indicates contamination from sewage because coliforms are a group of microorganisms fermenting sugars and producing carbondioxide. These microroganisms normally present human intestine and should not be present in dirnking water.